


Catheter Fixation Device

COMPOSITION:


Waterproof non-woven fabric, pressure sensitive adhesive, release paper, injection hook, fleece.

APPLICATION:

- 1. Indwelling needle fixation
- 2. Medical catheter and dressing fixation



Type A



Type B



Type I



Type J



Type E




Type C



Type H



Type F2



Type	Quantity/box	Quantity/case	Bag	Box size	Case size
All types	50pcs	900pcs	10x15cm	17x15.5x11cm	35x48x34.5cm

