

Gauze Non-woven Products




COMPOSITION:

Non-woven fabric/PU film,Acrylic adhesive,Cotton pad,Release paper,Protection paper

APPLICATION:

- 1.Protection for minor scratches, small area first and second degree burns, and minor wounds.
- 2.Prevention and treatment of bedsores and skin cracking.
- 3.Used for protection of donor site wounds and surgical incisions.
- 4.General protection for lacerations.



Gauze / Non-Woven	
Cutting Gauze	
12S-16S	21S-40S
30x20MESH	30x20MESH
10x10cm	10x10cm
20x20cm	20x20cm
30x30cm	30x30cm
40x40cm	40x40cm
50x50cm	50x50cm

Gauze / Non-Woven	
Sterile Gauze Swabs	
40S/26x18MESH	40S/30x20MESH
5cmx5cm	5cmx5cm
7.5cmx7.5cm	7.5cmx7.5cm
10x10cm-25P	10x10cm-25P
18x40cm	18x40cm
10x20cm	10x20cm
36x40cm	36x40cm