


Wound dressing

[Non-woven + PU]




COMPOSITION:

Non-woven fabric/PU film,Acrylic adhesive,
Cotton pad,Release paper,Protection paper

APPLICATION:

1. Protection for minor scratches, small area first and second degree burns, and minor wounds.
2. Prevention and treatment of bedsores and skin cracking.
3. Used for protection of donor site wounds and surgical incisions.



6x7cm



9.5x6.5cm




10x10cm



10x10cm




10x15cm



10x20cm



10x10cm



10x35cm

